



unique  **edge**[®]
MICROSURGICAL KNIVES

Performance, Precision *and* Protection


UniqueEdge[®] offers a new line of guarded microsurgical knives that are simple to use and make safety easy.

- Uses the same handle that you've always valued
- Requires no change in your surgical technique
- Retains maximal feel and absolute control
- Protects your staff from accidental injuries
- Protects the blade from accidental damage
- Part of your OSHA Exposure Control Plan (ECP)




To Use the Knife, Simply:


1 Twist



2 Retract



3 Lock



To Dispose of the Knife, Reverse Steps and Place in a Sharps Container


Guarded Stab Knives

ITEM NO.	WIDTH	DESCRIPTION
52-1501-G	15°	Full Handle
52-3001-G	30°	Full Handle



Guarded Crescent Knives

ITEM NO.	WIDTH	DESCRIPTION
54-1010-G	2.0mm	Angled, full handle, bevel up



Guarded Phaco Slit Knives

ITEM NO.	WIDTH	DESCRIPTION
52-2761-G	2.75mm	Angled, full handle, bevel up
52-2861-G	2.85mm	Angled, full handle, bevel up
52-2961-G	3.0mm	Angled, full handle, bevel up
52-3261-G	3.2mm	Angled, full handle, bevel up



Guarded Sideport Knives - MVR Design

ITEM NO.	WIDTH	DESCRIPTION
58-2010-G	1.0mm	Angled, full handle



Guarded MicroPhaco Slit Knives

ITEM NO.	WIDTH	DESCRIPTION
52-2263-G	2.2mm	Angled, full handle, double bevel



TEXTURED HANDLE

- Market proven design
- Enhanced instrument control
- Excellent tactile response

MATTE FINISHED BLADE

- Improves view of the operating field
- Eliminates reflective glare

PROTECTIVE GUARD*

- Protects from sharps injuries
- Fully retractable
- Secure integrated locking mechanism
- Simple, easy to use



To locate a UniqueEdge® Distributor or for ordering information, contact us at www.unique-tech-inc.com or by Phone at 610-775-9191

unique  **technologies inc.**

Mohnton PA USA Office: 610-775-9191 Fax: 610-775-9192

100% Made in the U.S.A.

ISO 9001-2000 certified & ISO 13485-2003 certified
CE Mark Authorized under MDD 93/42/EEC

* patent pending

MM102 R01

