



**UniBevel Knives –
THE STANDARD OF SHARPNESS FOR
THE MOST DEMANDING PROCEDURES**

- Jewel blade performance in a metal knife
- Super sharp for effortless penetration, clean incisions, and excellent control
- Ultra thin for minimal tissue damage
- Unique seamless, convex, bevel design
- Available in both single bevel Slit and double bevel Clear Corneal designs

**UniqueEdge Guarded Knives –
ABSOLUTE PRECISION AND PROTECTION**

- Same handle as our non-guarded knife for maximum comfort and control
- Easy to use sheath fully retracts for a completely unobstructed view
- Requires no change in surgical technique
- Helps eliminate accidental injuries to staff
- Protects the blade from damage
- Satisfies the OSHA Exposure Control Plan (ECP) and the European Council Directive 2010/32/EU

**To use a UniqueEdge Guarded Knife ,
simply:**



UniqueEdge knives are 100% made in the U.S.A., are satisfaction guaranteed, and available through a global network of authorized, professional distributors. To find a UniqueEdge distributor near you, please view our Distributor Locator, www.unique-tech-inc.com/distributorlocator.aspx

PRIVATE LABEL OPPORTUNITIES

We manufacture sterile and non-sterile microsurgical products for some of the largest and most trusted names in healthcare. We offer complete customization from product design to packaging. Contact us with your specific product needs.

WORLD CLASS QUALITY

Our quality management practices guide all of our manufacturing processes to ensure surgeons are using only the best performing products.

- Registered FDA Medical Device Manufacturer (US FDA Class 1 Exempt)
- ISO 9001-2008 certified
- ISO 13485-2012 certified
- ISO 13485-2003 certified

Our registration and certifications assure consistency and reliability of Unique's manufacturing processes and product quality. You can download these certificates from our website under the quality tab.



To find a UniqueEdge distributor, visit our Distributor Locator on our website, www.unique-tech-inc.com/distributorlocator.aspx

unique technologies inc.

Mohnton PA USA
Office: 610-775-9191 Fax: 610-775-9192

✦ www.unique-tech-inc.com



unique edge
MICROSURGICAL KNIVES

Precision you can see. Quality you can feel.

Precision you can see. Quality you can feel.

Unique Technologies, Inc. is a manufacturer and global distributor of precision microsurgical instruments. We develop and manufacture premium, single-use cutting instruments for some of the most trusted names in the microsurgical industry. Our UniqueEdge knives are known for being consistent, incredibly sharp, and the highest quality demanded by leading ophthalmic surgeons worldwide.

Clinicians know that minimizing risk and achieving positive patient outcomes are the most important issues to consider when choosing the right surgical instruments. UniqueEdge knives provide the quality, the precision, and the consistency required to achieve these goals.

PRECISION YOU CAN SEE

- Our vacuum-purified double melt, hardened and tempered stainless steel technology produces a uniform strength and ductility that is critical to creating a smooth, well-controlled incision
- Consistent blade dimensions
- Precise, flawless cutting edge
- Uniform, reproducible incisions
- Superior wound architecture



QUALITY YOU CAN FEEL

- Effortless tissue penetration
- Excellent tactile response
- Minimal drag
- Enhanced control

PERFORMANCE YOU CAN TRUST


- Accurate incisions
- Minimal wound leakage
- Tight self-sealing incisions
- Rapid healing
- Reduced potential for infection

TEXTURED HANDLE

- Surgeon proven design
- Comfortable secure grip
- Absolute instrument control

MATTE FINISHED BLADE

- Improves view of the operating field
- Eliminates reflective glare from operating room light sources
- Maximal visibility of the knife within the eye chamber



Slit Knives single bevel, bevel up

DESCRIPTION	ANGLED	STRAIGHT	GUARDED	THREADED TIP
1.4 mm	52-1463			
1.5 mm	52-1563			
1.8 mm	52-1863			
2.0 mm	52-2063		52-2063-G	
2.2 mm	52-2261		52-2261-G	
2.4 mm	52-2461		52-2461-G	
2.5 mm	52-2561		52-2561-G	
2.65 mm	52-2661		52-2661-G	52-2662
2.75 mm	52-2761	52-2701	52-2761-G	52-2762
2.8 mm	52-2808	52-2801	52-2808-G	
2.85 mm	52-2861		52-2861-G	52-2862
3.0 mm	52-2961	52-2901	52-2961-G	52-2962
3.2 mm	52-3261	52-3201	52-3261-G	52-3262
3.5 mm	52-3561		NA	
5.2 mm	52-5261		NA	



double bevel

DESCRIPTION	ANGLED	STRAIGHT	GUARDED
2.2 mm	52-2263		52-2263-G
2.4 mm	52-2431		
2.5 mm	52-2531		
2.65 mm	52-2631		
2.75 mm	52-2731		
2.8 mm	52-2830		
2.85 mm	52-2831		
3.0 mm	52-2931		
3.2 mm	52-3231		



UniBevel™ Slit Knives single bevel, bevel up

DESCRIPTION	ANGLED	STRAIGHT	GUARDED
1.0 mm	52-1040		52-1040-G
1.2 mm	52-1240		
1.4 mm	52-1440		
1.6 mm	52-1640		
1.8 mm	52-1840		52-1840-G
2.0 mm	52-2040		52-2040-G
2.2 mm	52-2240		52-2240-G
2.4 mm	52-2440		52-2440-G
2.5 mm	52-2540		52-2540-G
2.65 mm	52-2640		52-2640-G
2.75 mm	52-2740		52-2740-G
2.8 mm	52-2848		52-2848-G
2.85 mm	52-2840		
3.0 mm	52-3040		52-3040-G
3.2 mm	52-3240		52-3240-G



AK Slit Knives single bevel

DESCRIPTION	ANGLED	STRAIGHT	GUARDED	THREADED TIP
2.2 mm	52-2250			
2.4 mm	52-2450			
2.75 mm	52-2750			
2.8 mm	52-2851			
2.85 mm	52-2850			
3.0 mm	52-3050			
3.2 mm	52-3250			



double bevel

DESCRIPTION	ANGLED	STRAIGHT	GUARDED
2.4 mm	52-2452		
2.75 mm	52-2752		

UniqueEdge Knife and Blade Guide

All knives are full handle unless specified

We offer more than 200 UniqueEdge microsurgical knives and blades in a range of designs and sizes. We continuously add new knives and blades in response to the changing requirements of the surgical community and we welcome the opportunity to creatively address your new and demanding procedures. You can download our complete and up-to-date Product Catalog on our website at www.unique-tech-inc.com. For a description of each knife style, click on the "UniqueEdge" tab and the link for "Knives and Blades".

Stab Knives

DESCRIPTION	STRAIGHT	GUARDED	THREADED TIP
15°	52-1501	52-1501-G	52-1502
22.5°	52-2201	52-2201-G	52-2202
30°	52-3001	52-3001-G	52-3002
45°	52-4501	52-4501-G	



restricted depth

DESCRIPTION	STRAIGHT	GUARDED	THREADED TIP
15° 3 mm	52-1531		52-1532
15° 5 mm	52-1551		52-1552
22.5° 4 mm	52-2241		
30° 3.5 mm	52-3036		
30° 5 mm	52-3051		



MVR Sideport Knives

DESCRIPTION	ANGLED	STRAIGHT	GUARDED
0.8 mm (23G)	58-2008	58-2308	
1.0 mm	58-2010		58-2010-G
1.1 mm	58-2011		58-2011-G
1.2 mm (20G)	58-2012		58-2012-G
1.2 mm (20G)		58-2312	58-2312-G
1.4 mm (19G)	58-2014	58-2314	



Vitrectomy Knives

DESCRIPTION	ANGLED	STRAIGHT
19 gauge		51-1901
20 gauge	51-2061	51-2001
23 gauge	51-2361	51-2301
25 gauge	51-2531	51-2501



Crescent Knives

DESCRIPTION	ANGLED	STRAIGHT	GUARDED	THREADED TIP
2.0 mm		54-1000		54-5000
2.0 mm	54-1010		54-1010-G	54-5062
2.2 mm	54-1012		54-1012-G	
2.3 mm	54-1013			



single bevel, bevel down

DESCRIPTION	ANGLED	STRAIGHT	GUARDED
2.0 mm	54-1020		

double bevel

DESCRIPTION	ANGLED	STRAIGHT	GUARDED
1.25 mm	54-1060		
2.0 mm	54-1031		

Clear Corneal Knives

DESCRIPTION	ANGLED	GUARDED
2.5 mm	52-2563	
2.65 mm	52-2663	
2.75 mm	52-2763	52-2763-G
2.85 mm	52-2863	
3.0 mm	52-3063	
3.2 mm	52-3263	



double bevel, side cutting

DESCRIPTION	ANGLED	GUARDED
2.65 mm	52-2635	52-2635-G
2.75 mm	52-2735	
2.85 mm	52-2835	
3.0 mm	52-3035	
3.2 mm	52-3235	



UniBevel Clear Corneal Knives

DESCRIPTION	ANGLED	GUARDED
2.2 mm	52-2242	
2.4 mm	52-2442	
2.5 mm	52-2542	
2.65 mm	52-2642	
2.75 mm	52-2742	
2.85 mm	52-2842	
3.0 mm	52-3042	52-3042-G
3.2 mm	52-3242	52-3242-G



Disc/Spoon Knives

DESCRIPTION	ANGLED	GUARDED
2.2 mm	54-2261	
3.0 mm	54-3061	



Implant Knives

DESCRIPTION	ANGLED
3.2 mm	56-5361
3.5 mm	56-3561
4.0 mm	56-4061
5.2 mm	56-5261



Precision Depth Knives

DESCRIPTION	STRAIGHT	GUARDED
300 micron	52-0300	
350 micron	52-0350	
550 micron	52-0550	
600 micron	52-0600	



Trapezoid Knives

DESCRIPTION	ANGLED	GUARDED
1.2 x 1.4 mm	55-1214	
1.4 x 1.6 mm	55-1416	
1.5 x 2.0 mm	55-1520	
1.6 x 1.8 mm	55-1618	
1.8 x 2.2 mm	55-1822	
2.0 x 2.2 mm	55-2022	
2.5 x 2.9 mm	55-2522	
2.75 x 3.0 mm	55-2723	
2.75 x 3.2 mm	55-2724	
2.85 x 3.2 mm	55-2824	
2.85 x 3.5 mm	55-2825	
3.0 x 3.2 mm	55-2924	



Specialty Knives and Blades

DESCRIPTION	STRAIGHT	BLADE ONLY	THREADED TIPS
Needle Knife	58-6810		
#57 Scleral	58-5703	58-5700	58-5702
#64		58-6400	
#66 Lamellar	58-6603	58-6600	58-6602
#67		58-6700	
#69	58-6903	58-6900	58-6902
2.8 mm DB slit		52-2832F	
3.0 mm slit		52-2902F	
3.2 mm slit		52-3202F	
Myringotomy:			
Spear tip narrow shaft		58-7120	
Lance tip narrow shaft		58-7121	



We can guard almost any knife up to 3.2mm in size except for vitrectomy and implant knives. The guarded item numbers listed are for standard guarded knives that are carried in stock. New guarded knives may have minimum quantity requirements.

All knives are packaged 6 per box. If you don't find the knife you are looking for, please contact your UniqueEdge distributor or visit our Distributor Locator at: www.unique-tech-inc.com/distributorlocator.aspx

